



January 1, 2021

Dear Valued Customer,

In July of 2020, the Global Enteral Device Suppliers Association (GEDSA) announced the REVISED ENFit® Connector Conversion Schedule: U.S. and Legacy Connector Production Phase Out Dates (see below http://stayconnected.org/wp-content/uploads/2020/09/200914_REVISED_ENFit%C2%AE_Connector_Conversion_Schedule_U.S._Canada_Legacy_Connector_Production-Phase_Out_Dates-1-1.pdf). This schedule establishes a timeline for member manufactures to phase out legacy feeding devices and transitional adaptors to meet the requirements of the ISO standard 80369-3, commonly know as ENFit®, which goal is to minimize the risk of tubing misconnections.

| Date | Action | Example |
|--------------------------------|---|---|
| July 1 st , 2021 | Legacy feeding tubes and cross- application adaptors will no longer be manufactured. |  |
| January 1 st , 2022 | Transition sets and adaptors sold separately from other devices will no longer be manufactured. |  |

As an original member of GEDSA, and a worldwide leading provider of ENFit® products, GBUK Enteral is committed to this important patient safety initiative and the conversion timeframe. GBUK Enteral exclusively offers ENFit® products for your enteral needs.

To ensure that you and your institution are prepared for the transition to ENFit® please contact your local GBUK Enteral distributor or contact our Customer Service at (866)528-9422 or via our website; <https://www.gbukgroup.com/brands/gbuk-ental/>

Sincerely,

Kevin Maki
 Head of North American Sales
 GBUK Group Ltd