

ENFit® PHARMACY RESOURCE GUIDE

WHAT is ENFit? WHY was it created?

WHAT does ENFit mean for medication that goes through your feeding tube?

Feeding tube connectors have changed and are now designed to screw together to minimize the risks of tubing misconnections. These new connectors are called ENFit. Feeding tube users who have switched to ENFit should inform pharmacy and home tube feeding providers that ENFit syringes and medicine bottle adaptors are needed for medications that are delivered through your feeding tube.



FAQs:

- 1) Where can I obtain ENFit compatible supplies for my medication administration?

If you are having an ENFit tube placed, obtain the initial necessary supplies from your hospital before being discharged. If you have switched to ENFit and you receive your medication supplies from a local pharmacy or home care supplier, let them know you now have a feeding tube with an ENFit connection and will need ENFit syringes and medicine bottle adaptors. If they are unfamiliar with ENFit, you can share the Legacy to ENFit crosswalk and QR code, both shown below, so they can see what products to order. If you purchase your supplies from a medical supplier or online, you can use the Legacy to ENFit crosswalk to help identify the names of the ENFit compatible items to search for and the <https://stayconnected.org/pharmaguide/> webpage provides information on where to find those items.

- 2) Do ENFit supplies require a prescription?

Need for a prescription will vary depending on the supply item and even your insurance. If you are needing a prescription, contact your healthcare provider and let them know that you need a prescription that indicates ENFit on the order. You can share the Legacy to ENFit crosswalk below to help them know what the names of the needed items are. One of the benefits of the standardized ENFit connection is that ENFit products will connect with other ENFit products, no matter which brand. If you do have a particular brand preference, be sure to let your healthcare provider know so they can include the name on the prescription.

Legacy Enteral Syringe



ENFit Enteral Syringes



Legacy Press-In Bottle Adaptors



ENFit Press-In Bottle Adaptors



Legacy Universal / Stepped Bottle Adaptors



ENFit Universal / Stepped Bottle Adaptors



Legacy Pharmacy Bottle Cap Adaptors



ENFit Pharmacy Bottle Cap Adaptors



Legacy Medication Straws



ENFit Medication Straws



Other ENFit Pharmacy Supplies

ENFit Fill Cap Coupler



ENFit Syringe Caps



ENFit Transition Connector Being Used to Withdraw Medication



SHARE THIS CARD WITH YOUR PHARMACY

 **GEDSA**
stay connected[®]

For enteral medication administration with the ENFit system, please visit GEDSA's pharmacy guide webpage for information and product links

<https://stayconnected.org/pharmaguide>



SHARE THIS CARD WITH PATIENTS & CAREGIVERS

ENFit[®] Exchange
Facebook Group

To exchange ENFit information & products with other patients and caregivers, please join the ENFit Exchange Facebook Group



JOIN HERE