

February 7, 2022

**MOOG®**

Dear Valued Customer,

As you may be aware, on September 14<sup>th</sup>, 2020, the Global Enteral Device Supplier Association (GEDSA) released a revised conversion schedule indicating that member manufacturers will phase out legacy enteral feeding products and transition adapters. This transition will meet ISO standard 80369-3, commonly referred to as ENFit®.

As a member of GEDSA, Moog fully supports this essential safety initiative, and we are committed to removing the transitional connector from our products so that they are ENFit-only compatible. In an effort to better align with the transition schedules of other manufacturers in the enteral industry, beginning early in the second quarter of 2022 Moog will discontinue the manufacture of our feeding sets with transition connectors and will move forward with ENFit-only connectors on our Infinity® and Infinity® Orange® enteral feeding pump sets and gravity sets.

Because some inventory with transitional connectors is already in the field, you may not begin receiving the ENFit-only sets until the latter part of the second quarter of 2022 or shortly thereafter. In the interim, we strongly encourage you to prepare yourselves for this transition so that you will be able to meet the ENFit needs of your patients.

For a short period after Moog introduces the ENFit-only feeding sets, we will sell separately packaged transitional connectors to accommodate patients with non-ENFit feeding tubes as they are transitioned to ENFit feeding tubes. This offering of transitional connectors will be available for a limited time only. Other suppliers may also offer transition connectors for purchase. Check with your distribution channels regarding transition connector availability.

With the transition to ENFit-only enteral feeding sets, Moog will be changing current product code part numbers. Attachment A below includes a crosswalk of current ENFit product codes with transitional connectors to ENFit-only codes without transitional connectors. You will note that the only change to the product codes is the change from a “-A” (*Adapter*) suffix to a “-E” (*ENFit only*) suffix. As “-A” product inventories are depleted, we will fulfill orders with the new “-E” codes.

Moog is committed to the safety of our customers and patients. We believe the use of ENFit products will contribute to improved patient safety. If you have any questions regarding the ENFit initiative or how you can prepare for these changes, please contact your local Moog sales representative, customer service at (800) 970-2337, or email us at [clinicalsupport@moog.com](mailto:clinicalsupport@moog.com).

Sincerely,



John Bacon  
Sr. Product Manager  
Moog Medical

Infinity® and Infinity Orange® are trademarks of Zevex, Inc. ENFit® is a trademark of GEDSA

**Moog Medical**  
4314 Zevex Park Lane  
Salt Lake City, UT 84123 USA  
[www.moogmedical.com](http://www.moogmedical.com)

tel 801.264.1001  
toll free 800.970.2337  
fax 801.264.1051

## Existing Infinity Feeding Sets vs ENFit-Only Feeding Sets



Current Feeding set with transition connector

Future Feeding set without transition connector

## ENFit Feeding Set Crosswalk Codes

Current ENFit Administration Sets <i>with</i> Transitional Connectors		Future ENFit-Only Administration Sets <i>without</i> Transitional Connectors	
Product Description	Product Code	Product Description	Product Code
Infinity® Safety Screw Enteral Pump Delivery Set with ENFit® Connector System and Transitional Stepped Connector	INF0020-A	Infinity® Safety Screw Enteral Pump Delivery Set with ENFit® Connector	INF0020-E
Infinity® 500ml Enteral Pump Delivery Set with ENFit® Connector System and Transitional Stepped Connector	INF0500-A	Infinity® 500 mL Enteral Pump Delivery Set with ENFit® Connector	INF0500-E
Infinity® 1200ml Enteral Pump Delivery Set with ENFit® Connector System and Transitional Stepped Connector	INF1200-A	Infinity® 1200 mL Enteral Pump Delivery Set with ENFit® Connector	INF1200-E
Infinity® Orange® Screw Cap Enteral Pump Delivery Set with ENFit® Connector System and Transitional Stepped Connector	INFOCAP-A	Infinity® Orange® Screw Cap Enteral Pump Delivery Set with ENFit® Connector	INFOCAP-E
Infinity Orange 100mL Enteral Pump Delivery Set with ENFit® Connector System and Transitional Stepped Connector	INF0100-A	Infinity Orange 100mL Enteral Pump Delivery Set with ENFit® Connector	INF0100-E
1200ml Gravity Feeding Delivery Set with ENFit® Connector System and Transitional Stepped Connector	GR1200-A	1200mL Gravity Feeding Delivery Set with ENFit® Connector	GR1200-E